

## RCOR2 antibody

### Product Information

Catalog No.:	FNab07213
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May act as a component of a corepressor complex that represses transcription.

### Immunogen information

Immunogen:	REST corepressor 2
Synonyms:	REST corepressor 2 RCOR2
Observed MW:	58 kDa
Uniprot ID :	Q8IZ40

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab07213( RCOR2 Antibody) at dilution of 1:1000