

RECK antibody

Product Information

Catalog No.:	FNab07224
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Negatively regulates matrix metalloproteinase-9(MMP-9) by suppressing MMP-9 secretion and by direct inhibition of its enzymatic activity. RECK down-regulation by oncogenic signals may facilitate tumor invasion and metastasis. Appears to also regulate MMP-2 and MT1-MMP, which are involved in cancer progression.

Immunogen information

Immunogen:	reversion-inducing-cysteine-rich protein with kazal motifs
Synonyms:	Reversion-inducing cysteine-rich protein with Kazal motifs (hRECK) Suppressor of tumorigenicity 15 protein RECK ST15
Observed MW:	110 kDa
Uniprot ID :	O95980

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



293T cells were subjected to SDS PAGE followed by western blot with FNab07224(RECK antibody) at dilution of 1:1000