

RECQL4 antibody

Product Information

Catalog No.:	FNab07226
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


DNA-dependent ATPase. May modulate chromosome segregation.

Immunogen information


Immunogen:	RecQ protein-like 4
Synonyms:	ATP-dependent DNA helicase Q4 DNA helicase, RecQ-like type 4 (RecQ4) RTS RecQ protein-like 4 RECQL4 RECQ4
Observed MW:	145-150 kDa
Uniprot ID :	O94761

Application


Reactivity:	Human
Tested Application:	ELISA, WB, IF, IHC, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney tissue slide using FNab07226(RECQL4 Antibody) at dilution of 1:50



IP Result of anti-RECQL4 (IP:FNab07226, 4ug; Detection:FNab07226 1:1000) with HeLa cells lysate 2500ug.



HeLa cells were subjected to SDS PAGE followed by western blot with FNab07226(RECQL4 antibody) at dilution of 1:1000