

REG4 antibody

Product Information

Catalog No.:	FNab07233
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Calcium-independent lectin displaying mannose-binding specificity and able to maintain carbohydrate recognition activity in an acidic environment. May be involved in inflammatory and metaplastic responses of the gastrointestinal epithelium.

Immunogen information

Immunogen:	regenerating islet-derived family, member 4
Synonyms:	Regenerating islet-derived protein 4 (REG-4) Gastrointestinal secretory protein REG-like protein Regenerating islet-derived protein IV (Reg IV) REG4 GISP RELP
Observed MW:	18 kDa
Uniprot ID :	Q9BYZ8

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



SGC-7901 cells were subjected to SDS PAGE followed by western blot with FNab07233(REG4 antibody) at dilution of 1:300