

RFTN1 antibody

Product Information

Catalog No.:	FNab07255
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May play a pivotal role in the formation and/or maintenance of lipid rafts. May regulate B-cell antigen receptor-mediated signaling.

Immunogen information


Immunogen:	raftlin, lipid raft linker 1
Synonyms:	Raftlin Cell migration-inducing gene 2 protein Raft-linking protein RFTN1 KIAA0084
Observed MW:	63-70 kDa
Uniprot ID :	Q14699

Application

Reactivity:	Human
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human tonsillitis slide using FNab07255(RFTN1 Antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab07255(RFTN1 Antibody) at dilution of 1:2000