

RGS11 antibody

Product Information

Catalog No.:	FNab07265
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Inhibits signal transduction by increasing the GTPase activity of G protein alpha subunits thereby driving them into their inactive GDP-bound form.

Immunogen information

Immunogen:	regulator of G-protein signaling 11
Synonyms:	Regulator of G-protein signaling 11 (RGS11) RGS11
Observed MW:	53 kDa
Uniprot ID :	O94810

Application

Reactivity:	Human, Rat, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



HL-60 cells were subjected to SDS PAGE followed by western blot with FNab07265(RGS11 Antibody) at dilution of 1:1000