

RHBDL2 antibody

Product Information

Catalog No.:	FNab07278
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Involved in regulated intramembrane proteolysis and the subsequent release of functional polypeptides from their membrane anchors. Known substrate: EFNB3.

Immunogen information

Immunogen:	rhomboid, veinlet-like 2
Synonyms:	Rhomboid-related protein 2 (RRP2) Rhomboid-like protein 2 Rhomboid-related protein 2, N-terminal fragment (NTF) Rhomboid-related protein 2, C-terminal fragment (CTF) RHBDL2
Observed MW:	35 kDa
Uniprot ID :	Q9NX52

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:1000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab07278(RHBDL2 antibody) at dilution of 1:500