

RHCG antibody

Product Information

Catalog No.:	FNab07279
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Functions as an electroneutral and bidirectional ammonium transporter. May regulate transepithelial ammonia secretion.

Immunogen information


Immunogen:	Rh family, C glycoprotein
Synonyms:	Ammonium transporter Rh type C Rh glycoprotein kidney Rhesus blood group family type C glycoprotein (Rh family type C glycoprotein, Rh type C glycoprotein) Tumor-related protein DRC2 RHCG C15orf6 CDRC2 PDRC2 RHGK
Observed MW:	52-55 kDa
Uniprot ID :	Q9UBD6

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas cancer using FNab07279(RHCG antibody) at dilution of 1:50



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab07279(RHCG antibody) at dilution of 1:300