

## RHOBTB3 antibody

### Product Information

Catalog No.:	FNab07283
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Rab9-regulated ATPase required for endosome to Golgi transport. Involved in transport vesicle docking at the Golgi complex, possibly by participating in release M6PRBP1/TIP47 from vesicles to permit their efficient docking and fusion at the Golgi. Specifically binds Rab9, but not other Rab proteins. Has low intrinsic ATPase activity due to autoinhibition, which is relieved by Rab9.

### Immunogen information


Immunogen:	Rho-related BTB domain containing 3
Synonyms:	Rho-related BTB domain-containing protein 3 RHOBTB3 KIAA0878
Observed MW:	69 kDa
Uniprot ID :	O94955

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human placenta tissue slide using FNab07283(RHOBTB3 Antibody) at dilution of 1:50



K-562 cells were subjected to SDS PAGE followed by western blot with FNab07283(RHOBTB3 antibody) at dilution of 1:1000