

## RHOC antibody

### Product Information

Catalog No.:	FNab07284
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Regulates a signal transduction pathway linking plasma membrane receptors to the assembly of focal adhesions and actin stress fibers. Serves as a microtubule-dependent signal that is required for the myosin contractile ring formation during cell cycle cytokinesis. Regulates apical junction formation in bronchial epithelial cells.

### Immunogen information

Immunogen:	ras homolog gene family, member C
Synonyms:	Rho-related GTP-binding protein RhoC Rho cDNA clone 9 (h9) RHOC ARH9 ARHC
Observed MW:	22 kDa
Uniprot ID :	P08134

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:200-1:2000
Image:	



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with FNab07284(RHOC antibody) at dilution of 1:500