

RHOF antibody

Product Information

Catalog No.:	FNab07285
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Plasma membrane-associated small GTPase which cycles between an active GTP-bound and an inactive GDP-bound state. Causes the formation of thin, actin-rich surface projections called filopodia. Functions cooperatively with CDC42 and Rac to generate additional structures, increasing the diversity of actin-based morphology.

Immunogen information


Immunogen:	ras homolog gene family, member F(in filopodia)
Synonyms:	Rho-related GTP-binding protein RhoF Rho family GTPase Rif Rho in filopodia RHOF ARHF RIF
Observed MW:	24 kDa, 19 kDa
Uniprot ID :	Q9HBH0

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer using FNa07285(RHOF antibody) at dilution of 1:50



COLO 320 cells were subjected to SDS PAGE followed by western blot with FNa07285(RHOF antibody) at dilution of 1:500