

RHPN1 antibody

Product Information

Catalog No.:	FNab07289
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Has no enzymatic activity. May serve as a target for Rho, and interact with some cytoskeletal component upon Rho binding or relay a Rho signal to other molecules(By similarity).

Immunogen information

Immunogen:	rhophilin, Rho GTPase binding protein 1
Synonyms:	Rhophilin-1 GTP-Rho-binding protein 1 RHPN1 KIAA1929
Observed MW:	71 kDa
Uniprot ID :	Q8TCX5

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



human brain tissue were subjected to SDS PAGE followed by western blot with FNab07289(RHPN1 antibody) at dilution of 1:400