

## RIMS1 antibody

### Product Information

Catalog No.:	FNab07300
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background

None

### Immunogen information

Immunogen:	regulating synaptic membrane exocytosis 1
Synonyms:	Regulating synaptic membrane exocytosis protein 1 Rab-3-interacting molecule 1 (RIM 1) Rab-3-interacting protein 2 RIMS1 KIAA0340 RAB3IP2 RIM1
Observed MW:	189 kDa
Uniprot ID :	Q86UR5

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:1000-1:4000
Image:	

mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab07300(RIM1 antibody) at dilution of 1:2000

