

## RIMS2-Specific antibody

### Product Information

Catalog No.:	FNab07301
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	regulating synaptic membrane exocytosis 2
Synonyms:	Regulating synaptic membrane exocytosis protein 2 Rab-3-interacting molecule 2 (RIM 2) Rab-3-interacting protein 3 RIMS2 KIAA0751 RAB3IP3 RIM2
Observed MW:	133 kDa
Uniprot ID :	Q9UQ26

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



L02 cells were subjected to SDS PAGE followed by western blot with FNab07301(RIM2-Specific antibody) at dilution of 1:400