

RIN3 antibody

Product Information

Catalog No.:	FNab07307
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Ras effector protein that functions as a guanine nucleotide exchange(GEF) for RAB5B and RAB31, by exchanging bound GDP for free GTP. Required for normal RAB31 function.

Immunogen information

Immunogen:	Ras and Rab interactor 3
Synonyms:	Ras and Rab interactor 3 Ras interaction/interference protein 3 RIN3
Observed MW:	108 kDa
Uniprot ID :	Q8TB24

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:200-1:1000
Image:	



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab07307(RIN3 antibody) at dilution of 1:300