

RIOK1 antibody

Product Information

Catalog No.:	FNab07310
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	RIO kinase 1
Synonyms:	Serine/threonine-protein kinase RIO1 RIO kinase 1 RIOK1 RIO1
Observed MW:	75-80 kDa
Uniprot ID :	Q9BRS2

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Hela cells were subjected to SDS PAGE followed by western blot with FNab07310(RIOK1 antibody) at dilution of 1:1000