

CLVS2 antibody

Product Information

Catalog No.:	FNab07319
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	retinaldehyde binding protein 1-like 2
Synonyms:	Clavesin-2 Retinaldehyde-binding protein 1-like 2 clathrin vesicle-associated Sec14 protein 2 CLVS2 C6orf212 C6orf213 RLBP1L2
Observed MW:	38 kDa
Uniprot ID :	Q5SYC1

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab07319(RLBP1L2 antibody) at dilution of 1:1500