

RNASE13 antibody

Product Information

Catalog No.:	FNab07326
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)


Background

Immunogen information


Immunogen:	ribonuclease, RNase A family, 13 (non-active)
Synonyms:	Probable inactive ribonuclease-like protein 13 RNASE13
Observed MW:	18 kDa
Uniprot ID :	Q5GAN3

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:100 - 1:200; IF: 1:10 - 1:100
Image:	



Immunohistochemistry of paraffin-embedded human stomach tissue slide using FNab07326(RNASE13 Antibody) at dilution of 1:100



A431 cells were subjected to SDS PAGE followed by western blot with FNab07326(RNASE13 antibody) at dilution of 1:1000