

## RNF112 antibody

### Product Information

Catalog No.:	FNab07336
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	ring finger protein 112
Synonyms:	RING finger protein 112 Brain finger protein Zinc finger protein 179 RNF112 BFP ZNF179
Observed MW:	68 kDa
Uniprot ID :	Q9ULX5

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



mouse kidney tissue were subjected to SDS PAGE followed by western blot with FNab07336( RNF112 Antibody) at dilution of 1:600