

MGRN1 antibody

Product Information

Catalog No.:	FNab07342
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	mahogunin, ring finger 1
Synonyms:	E3 ubiquitin-protein ligase MGRN1 Mahogunin RING finger protein 1 RING finger protein 156 RING-type E3 ubiquitin transferase MGRN1 MGRN1 KIAA0544 RNF156
Observed MW:	58-64 kDa
Uniprot ID :	O60291

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lung cancer using FNab07342(MGRN1 antibody) at dilution of 1:100



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab07342(MGRN1 antibody) at dilution of 1:500