

RNF25 antibody

Product Information

Catalog No.:	FNab07353
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


E3 ubiquitin-protein ligase that mediates ubiquitination and subsequent proteasomal degradation of NKD2(By similarity). Stimulates transcription mediated by NF-kappa-B.

Immunogen information


Immunogen:	ring finger protein 25
Synonyms:	E3 ubiquitin-protein ligase RNF25 RING finger protein 25 RING finger protein AO7 RNF25
Observed MW:	51 kDa
Uniprot ID :	Q96BH1

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human cervical cancer tissue slide using FNab07353(RNF25 Antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab07353(RNF25 Antibody) at dilution of 1:3000