

## MRM3 antibody

### Product Information

Catalog No.:	FNab07368
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	RNA methyltransferase like 1
Synonyms:	rRNA methyltransferase 3, mitochondrial 16S rRNA (guanosine(1370)-2'-O)-methyltransferase 16S rRNA [Gm1370] 2'-O-methyltransferase RNA methyltransferase-like protein 1 MRM3 RNMTL1
Observed MW:	47 kDa
Uniprot ID :	Q9HC36

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver using FNab07368(RNMTL1 antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with FNab07368(RNMTL1 antibody) at dilution of 1:500