

## ROR1 antibody

### Product Information

Catalog No.:	FNab07386
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Has very low kinase activity in vitro and is unlikely to function as a tyrosine kinase in vivo. May act as a receptor for wnt ligand WNT5A which may result in the inhibition of WNT3A-mediated signaling.

### Immunogen information


Immunogen:	receptor tyrosine kinase-like orphan receptor 1
Synonyms:	Inactive tyrosine-protein kinase transmembrane receptor ROR1 Neurotrophic tyrosine kinase, receptor-related 1 ROR1 NTRKR1
Observed MW:	130-140 kDa
Uniprot ID :	Q01973

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab07386(ROR1 Antibody) at dilution of 1:50



K-562 cells were subjected to SDS PAGE followed by western blot with FNab07386(ROR1 antibody) at dilution of 1:500