

RORC antibody

Product Information


Catalog No.:	FNab07389
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	RAR-related orphan receptor C
Synonyms:	Nuclear receptor ROR-gamma Nuclear receptor RZR-gamma Nuclear receptor subfamily 1 group F member 3 RAR-related orphan receptor C Retinoid-related orphan receptor-gamma RORC NR1F3 RORG RZRG
Observed MW:	58 kDa
Uniprot ID :	P51449

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500 - 1:2000
Image:	



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab07389(RORC antibody) at dilution of 1:1000