

RPA3 antibody

Product Information

Catalog No.:	FNab07393
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	replication protein A3, 14kDa
Synonyms:	Replication protein A 14 kDa subunit (RP-A p14) Replication factor A protein 3 (RF-A protein 3) RPA3 REPA3 RPA14
Observed MW:	None
Uniprot ID :	P35244

Application

Reactivity:	Human
Tested Application:	ELISA
Recommended dilution:	None