

RPAP3 antibody

Product Information

Catalog No.:	FNab07398
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Forms an interface between the RNA polymerase II enzyme and chaperone/scaffolding protein, suggesting that it is required to connect RNA polymerase II to regulators of protein complex formation.

Immunogen information

Immunogen:	RNA polymerase II associated protein 3
Synonyms:	RNA polymerase II-associated protein 3 RPAP3
Observed MW:	76 kDa
Uniprot ID :	Q9H6T3

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab07398(RPAP3 antibody) at dilution of 1:1000