

## RPGRIP1 antibody

### Product Information

Catalog No.:	FNab07403
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May function as scaffolding protein. Required for normal location of RPGR at the connecting cilium of photoreceptor cells. Required for normal disk morphogenesis and disk organization in the outer segment of photoreceptor cells and for survival of photoreceptor cells.

### Immunogen information

Immunogen:	retinitis pigmentosa GTPase regulator interacting protein 1
Synonyms:	X-linked retinitis pigmentosa GTPase regulator-interacting protein 1 (RPGR-interacting protein 1) RPGRIP1
Observed MW:	70 kDa, 140-150 kDa
Uniprot ID :	Q96KN7

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab07403(RPGRIP1 antibody) at dilution of 1:400