

RPL12 antibody

Product Information

Catalog No.:	FNab07411
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Binds directly to 26S ribosomal RNA.

Immunogen information

Immunogen:	ribosomal protein L12
Synonyms:	Large ribosomal subunit protein uL11 60S ribosomal protein L12 RPL12
Observed MW:	21 kDa
Uniprot ID :	P30050

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab07411(RPL12 antibody) at dilution of 1:400