

RPL18A antibody

Product Information

Catalog No.:	FNab07418
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	ribosomal protein L18a
Synonyms:	Large ribosomal subunit protein eL20 60S ribosomal protein L18a RPL18A
Observed MW:	21 kDa
Uniprot ID :	Q02543

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:3000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human brain using FNab07418(RPL18A antibody) at dilution of 1:100



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab07418(RPL18A antibody) at dilution of 1:500