

RPL22 antibody

Product Information

Catalog No.:	FNab07421
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	ribosomal protein L22
Synonyms:	Large ribosomal subunit protein eL22 60S ribosomal protein L22 EBER-associated protein (EAP) Epstein-Barr virus small RNA-associated protein Heparin-binding protein HBp15 RPL22
Observed MW:	15-18 kDa
Uniprot ID :	P35268

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



A431 cells were subjected to SDS PAGE followed by western blot with FNab07421(RPL22 Antibody) at dilution of 1:600