

## RPL27A antibody

### Product Information

Catalog No.:	FNab07427
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	ribosomal protein L27a
Synonyms:	Large ribosomal subunit protein uL15 60S ribosomal protein L27a RPL27A
Observed MW:	18 kDa
Uniprot ID :	P46776

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



A549 cells were subjected to SDS PAGE followed by western blot with FNab07427( RPL27A Antibody) at dilution of 1:300