

## RPL3 antibody

### Product Information

Catalog No.:	FNab07429
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


The L3 protein is a component of the large subunit of cytoplasmic ribosomes.

### Immunogen information


Immunogen:	ribosomal protein L3
Synonyms:	Large ribosomal subunit protein uL3 60S ribosomal protein L3 HIV-1 TAR RNA-binding protein B (TARBP-B) RPL3
Observed MW:	46 kDa
Uniprot ID :	P39023

### Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, IHC, IF, IP, WB
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human prostate cancer using FNab07429(RPL3 antibody) at dilution of 1:10



IP Result of anti-RPL3 (IP:FNab07429, 4ug; Detection:FNab07429 1:500) with mouse kidney tissue lysate 4000ug.



mouse placenta tissue were subjected to SDS PAGE followed by western blot with FNab07429(RPL3 antibody) at dilution of 1:1000