

RPP25 antibody

Product Information

Catalog No.:	FNab07453
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Component of ribonuclease P, a protein complex that generates mature tRNA molecules by cleaving their 5'-ends. Also a component of RNase MRP. This subunit binds to RNA.

Immunogen information


Immunogen:	ribonuclease P/MRP 25kDa subunit
Synonyms:	Ribonuclease P protein subunit p25 (RNase P protein subunit p25) RPP25
Observed MW:	25 kDa
Uniprot ID :	Q9BUL9

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human cervical cancer tissue slide using FNab07453(RPP25 Antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab07453(RPP25 antibody) at dilution of 1:300