

## RPS16 antibody

### Product Information

|               |                                                                                                                |
|---------------|----------------------------------------------------------------------------------------------------------------|
| Catalog No.:  | FNab07462                                                                                                      |
| Size:         | 100µg                                                                                                          |
| Form:         | liquid                                                                                                         |
| Purification: | Immunogen affinity purified                                                                                    |
| Purity:       | ≥95% as determined by SDS-PAGE                                                                                 |
| Host:         | Rabbit                                                                                                         |
| Clonality:    | polyclonal                                                                                                     |
| Clone ID:     | None                                                                                                           |
| IsoType:      | IgG                                                                                                            |
| Storage:      | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

### Background


None

### Immunogen information


|              |                                                                     |
|--------------|---------------------------------------------------------------------|
| Immunogen:   | ribosomal protein S16                                               |
| Synonyms:    | Small ribosomal subunit protein uS9 40S ribosomal protein S16 RPS16 |
| Observed MW: | 16 kDa                                                              |
| Uniprot ID : | P62249                                                              |

### Application

|                       |                                                   |
|-----------------------|---------------------------------------------------|
| Reactivity:           | Human, Mouse, Rat                                 |
| Tested Application:   | ELISA, WB, IHC, IF                                |
| Recommended dilution: | WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200 |
| Image:                |                                                   |



Immunohistochemistry of paraffin-embedded human breast cancer using FNab07462(RPS16 antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab07462(RPS16 antibody) at dilution of 1:500