

RPS27L antibody

Product Information

Catalog No.:	FNab07472
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	ribosomal protein S27-like
Synonyms:	Ribosomal protein eS27-like 40S ribosomal protein S27-like Small ribosomal subunit protein eS27-like RPS27L
Observed MW:	10 kDa
Uniprot ID :	Q71UM5

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



SH-SY5Y cells were subjected to SDS PAGE
followed by western blot with
FNab07472(RPS27L antibody) at dilution of 1:300