

RPS5 antibody

Product Information

Catalog No.:	FNab07481
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	ribosomal protein S5
Synonyms:	Small ribosomal subunit protein uS7 40S ribosomal protein S5 Small ribosomal subunit protein uS7, N-terminally processed Alternative names: 40S ribosomal protein S5, N-terminally processed RPS5
Observed MW:	23 kDa
Uniprot ID :	P46782

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab07481(RPS5 antibody) at dilution of 1:600