

RPS9 antibody

Product Information

Catalog No.:	FNab07486
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	ribosomal protein S9
Synonyms:	Small ribosomal subunit protein uS4 40S ribosomal protein S9 RPS9
Observed MW:	25 kDa
Uniprot ID :	P46781

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney using FNab07486(RPS9 antibody) at dilution of 1:100



mouse uterus tissue were subjected to SDS PAGE followed by western blot with FNab07486(RPS9 antibody) at dilution of 1:500