

RTCA antibody

Product Information

Catalog No.:	FNab07516
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	RNA terminal phosphate cyclase domain 1
Synonyms:	RNA 3'-terminal phosphate cyclase (RNA cyclase, RNA-3'-phosphate cyclase) RNA terminal phosphate cyclase domain-containing protein 1 (RTC domain-containing protein 1) RTCA RPC RPC1 RTC1 RTCD1
Observed MW:	41-43 kDa
Uniprot ID :	O00442

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



K-562 cells were subjected to SDS PAGE followed by western blot with FNab07516(RTCD1 antibody) at dilution of 1:500