

## RTN1-C antibody

### Product Information

Catalog No.:	FNab07521
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Human RTN1 (also known as NSP) gene is expressed predominantly in neuroendocrine tissues, and produces three isoforms termed RTN1-A (NSP-A), RTN1-B (NSP-B), and RTN1-C (NSP-C). This antibody is raised against the full-length protein of human RTN1-C.

### Immunogen information


Immunogen:	reticulon 1
Synonyms:	Reticulon-1 Neuroendocrine-specific protein RTN1 NSP
Observed MW:	23 kDa
Uniprot ID :	Q16799

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human brain using FNab07521(RTN1 antibody) at dilution of 1:100



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab07521(RTN1 antibody) at dilution of 1:800