

RTN2 antibody

Product Information

Catalog No.:	FNab07522
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	reticulon 2
Synonyms:	Reticulon-2 Neuroendocrine-specific protein-like 1 (NSP-like protein 1) Neuroendocrine-specific protein-like I (NSP-like protein I, NSPLI) RTN2 NSPL1
Observed MW:	20-23 kDa
Uniprot ID :	O75298

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human skeletal muscle slide using FNab07522(RTN2 Antibody) at dilution of 1:50



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab07522(RTN2 Antibody) at dilution of 1:1000