

## RTN4 antibody

### Product Information

Catalog No.:	FNab07524
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	reticulon 4
Synonyms:	Reticulon-4 Foocen Neurite outgrowth inhibitor (Nogo protein) Neuroendocrine-specific protein (NSP) Neuroendocrine-specific protein C homolog RTN-x Reticulon-5 RTN4 KIAA0886 NOGO
Observed MW:	45 kDa, 190-210 kDa
Uniprot ID :	Q9NQC3

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, FC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab07524(RTN4 Antibody) at dilution of 1:50



A375 cells were subjected to SDS PAGE followed by western blot with FNab07524(RTN4 antibody) at dilution of 1:1000