

## RUBCN antibody

### Product Information

Catalog No.:	FNab07527
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	KIAA0226
Synonyms:	Run domain Beclin-1-interacting and cysteine-rich domain-containing protein (Rubicon) Beclin-1 associated RUN domain containing protein (Baron) RUBCN KIAA0226
Observed MW:	110-130 kDa
Uniprot ID :	Q92622

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis using FNab07527(KIAA0226 antibody) at dilution of 1:50



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab07527(KIAA0226 antibody) at dilution of 1:1000