

RUFY2 antibody

Product Information

Catalog No.:	FNab07529
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)


Background

Immunogen information


Immunogen:	RUN and FYVE domain containing 2
Synonyms:	RUN and FYVE domain-containing protein 2 Rab4-interacting protein related RUFY2 KIAA1537 RABIP4R
Observed MW:	75 kDa
Uniprot ID :	Q8WXA3

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded human brain cancer using FNab07529(RUFY2 Antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with FNab07529(RUFY2 antibody) at dilution of 1:1000