

## RUNDC3A antibody

### Product Information

Catalog No.:	FNab07531
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May act as an effector of RAP2A in neuronal cells.

### Immunogen information


Immunogen:	RUN domain containing 3A
Synonyms:	RUN domain-containing protein 3A Rap2-interacting protein 8 (RPIP-8) RUNDC3A RAP2IP RPIP8
Observed MW:	66 kDa
Uniprot ID :	Q59EK9

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis using FNab07531(RUNDC3A antibody) at dilution of 1:100



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab07531(RUNDC3A antibody) at dilution of 1:2000