

RUNDC3B antibody

Product Information

Catalog No.:	FNab07532
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	RUN domain containing 3B
Synonyms:	RUN domain-containing protein 3B Rap2-binding protein 9 Rap2-interacting protein 9 (RPIP-9) RUNDC3B RPIB9 RPIP9
Observed MW:	45 kDa
Uniprot ID :	Q96NL0

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab07532(RUNDC3B antibody) at dilution of 1:500