

RYR2 antibody

Product Information

Catalog No.:	FNab07546
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	ryanodine receptor 2(cardiac)
Synonyms:	Ryanodine receptor 2 (RYR-2, RyR2, hRYR-2) Cardiac muscle ryanodine receptor Cardiac muscle ryanodine receptor-calcium release channel Type 2 ryanodine receptor RYR2
Observed MW:	250 kDa
Uniprot ID :	Q92736

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:100-1:400
Image:	



Immunohistochemistry of paraffin-embedded human heart tissue slide using FNab07546(RYR2 Antibody) at dilution of 1:200



HeLa cells were subjected to SDS PAGE followed by western blot with FNab07546(RYR2 antibody) at dilution of 1:500