

## S100A3 antibody

### Product Information

Catalog No.:	FNab07559
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Binds both calcium and zinc. May be involved in calcium-dependent cuticle cell differentiation, hair shaft and hair cuticular barrier formation.

### Immunogen information

Immunogen:	S100 calcium binding protein A3
Synonyms:	Protein S100-A3 Protein S-100E S100 calcium-binding protein A3 S100A3 S100E
Observed MW:	12 kDa
Uniprot ID :	P33764

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:200-1:1000
Image:	



Hela cells were subjected to SDS PAGE followed by western blot with FNab07559(S100A3 antibody) at dilution of 1:500