

## S100A9 antibody

### Product Information


Catalog No.:	FNab07565
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information


Immunogen:	S100 calcium binding protein A9
Synonyms:	Protein S100-A9 Calgranulin-B Calprotectin L1H subunit Leukocyte L1 complex heavy chain Migration inhibitory factor-related protein 14 (MRP-14, p14) S100 calcium-binding protein A9 S100A9 CAGB CFAG MRP14
Observed MW:	14 kDa
Uniprot ID :	P06702

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab07565(S1009 Antibody) at dilution of 1:200



human heart tissue were subjected to SDS PAGE followed by western blot with FNab07565(S100A9 antibody) at dilution of 1:500