

## SAA1 antibody

### Product Information

Catalog No.:	FNab07572
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Major acute phase protein.

### Immunogen information


Immunogen:	serum amyloid A1
Synonyms:	Serum amyloid A-1 protein (SAA) Amyloid protein A Alternative names: Amyloid fibril protein AA Serum amyloid protein A(2-104) Serum amyloid protein A(3-104) Serum amyloid protein A(2-103) Serum amyloid protein A(2-102) Serum amyloid protein A(4-101) SAA1
Observed MW:	14 kDa
Uniprot ID :	P0DJ18

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	



Immunohistochemical of paraffin-embedded human liver using FNaB07572(SAA1 antibody) at dilution of 1:50



Rat liver tissue were subjected to SDS PAGE followed by western blot with FNaB07572(SAA1 antibody) at dilution of 1:1000