

SAA4 antibody

Product Information

Catalog No.:	FNab07574
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Major acute phase reactant. Apolipoprotein of the HDL complex.

Immunogen information


Immunogen:	serum amyloid A4, constitutive
Synonyms:	Serum amyloid A-4 protein Constitutively expressed serum amyloid A protein (C-SAA) SAA4 CSAA
Observed MW:	14 kDa
Uniprot ID :	P35542

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver cancer using FNab07574(SAA4 antibody) at dilution of 1:50



Hela cells were subjected to SDS PAGE followed by western blot with FNab07574(SAA4 Antibody) at dilution of 1:1000